



OneTouch Reveal® web app

Helps you get more meaning out of your results

- Tracks your blood glucose, food, insulin dosing and activity data over time
- Suggests personalized tips to help you stay on track and in range*
- Reports can be printed, ready for discussion at your next appointment with your healthcare professional
- Can be used on PC and Mac computers at OneTouchReveal.com

*Be sure to talk to your healthcare professional about the actions that are right for you.

How to setup your OneTouch Reveal® web app

Before you begin, confirm that cookies are enabled in your browser.



Wireless Option

1. Connect wireless devices and activate your account

- Install the OneTouch Reveal® mobile app on your mobile device
- Open the mobile app and follow the on-screen instructions to setup your account and pair your OneTouch Verio Flex™ meter
- Check your email inbox for a link to confirm your email address

For more details, refer to the OneTouch Reveal® mobile app [User Guide](#).

2. Sign in to the OneTouch Reveal® web app on your computer

- Go to [OneTouchReveal.com](#) and enter the same **Username** and **Password** that you used in the OneTouch Reveal® mobile app
- Click **Sign In** and follow the on-screen instructions

OR

Cable Option

1. Register and activate your account

- Go to [OneTouchReveal.com](#), click on **Create a New Account** and follow the prompts
- Check your email inbox for a link to confirm your email address

2. Install Data Transfer Tool

- Unplug all OneTouch® USB cables from your computer
- **Sign In** to your OneTouch Reveal® account and click **help** from the top right of any screen
- Click the **Download** button under **Data Transfer Tool** and follow the prompts
- Select the appropriate devices
- Click **Finish**. After installation, you may be prompted to re-start your computer.

3. Transfer data from your OneTouch® meter

- Connect your meter to the computer with your OneTouch® USB cable (*cable sold separately*).
- Open the OneTouch Reveal® Data Transfer Tool and select **Detect Now**. When the device is found, a confirmation screen will appear.
- **Sign In** to the OneTouch Reveal® web application
- Click **Confirm** to transfer data. When the transfer completes, click **OK**

You're all set! You can now access your data with the OneTouch Reveal® web app at [OneTouchReveal.com](#)

ONE TOUCH®
Every Touch is a Step Forward™



View your blood glucose readings

- Click the **Logbook** tab to see your blood glucose readings
- Click the **Patterns** tab to view a 14-day summary

Patient Summary: 5/5/2016 - 5/18/2016

Name: Patient 3, Kate Patient ID: TTCP3 Overall Target: 70 - 180 mg/dL
 DOB: 5/1/1987 Clinic: Trident Training Clinic Before Meal Target: 70 - 130 mg/dL
 Gender: Female Clinician: Chris Mack, MD, TIC After Meal Target: 70 - 180 mg/dL

Warning: Reports are intended for use by healthcare professionals with experience treating patients with diabetes only.

Alerts:

- Low Alerts: Disabled, 16 glucose readings were below target (12.0%)
- High Alerts: Enabled, 14 glucose readings were above target (10.0%)
- Time Alerts: Disabled, 14 glucose readings were above target (10.0%)
- Time Alerts: Disabled, 14 glucose readings were above target (10.0%)

Comprehensive Statistics

	Current Date Range	Previous Date Range	Change %
Mean Glucose (mmol/L)	6.86	6.86	0%
Mean Glucose (mg/dL)	123.3	123.3	0%
Standard Deviation	1.84	1.84	0%
Max	14.4	14.4	0%
Min	3.3	3.3	0%
Time in Range	54%	54%	0%

Comprehensive Statistics

	Before Target	After Target
Low Alerts	0	0
High Alerts	0	0
Time Alerts	0	0
Time Alerts	0	0

Download by Time of Day

Download by Time of Day

Logbook: 5/17/2016 - 5/30/2016

Name: Patient 3, Kate Patient ID: TTCP3 Overall Target: 70 - 180 mg/dL
 DOB: 5/1/1987 Clinic: Trident Training Clinic Before Meal Target: 70 - 130 mg/dL
 Gender: Female Clinician: Chris Mack, MD, TIC After Meal Target: 70 - 180 mg/dL

Warning: Reports are intended for use by healthcare professionals with experience treating patients with diabetes only.

Date	Before Breakfast	After Breakfast	Before Lunch	After Lunch	Before Dinner	After Dinner	Bedtime
5/17/2016	4.6 @ 7:30 am	115 @ 7:30 am	83 @ 11:30 am	111 @ 12:30 pm	109 @ 7:50 pm	78 @ 10:50 pm	77 @ 10:20 pm
5/18/2016	4.6 @ 7:30 am	91 @ 7:00 am	88 @ 9:30 am	94 @ 12:30 pm	105 @ 6:30 pm	84 @ 8:47 pm	100 @ 10:00 pm
5/19/2016	2.6 @ 12:00 am	100 @ 8:40 am	112 @ 7:24 am	114 @ 1:50 pm	180 @ 8:20 pm	146 @ 10:00 pm	85 @ 10:00 pm
5/20/2016	4.6 @ 7:30 am	140 @ 4:42 am	112 @ 7:24 am	114 @ 1:50 pm	180 @ 8:20 pm	146 @ 10:00 pm	85 @ 10:00 pm
5/21/2016	4.6 @ 7:30 am	179 @ 6:27 am	88 @ 9:10 am	87 @ 12:07 pm	191 @ 8:13 pm	144 @ 8:28 pm	176 @ 10:50 pm
5/22/2016	4.6 @ 7:30 am	83 @ 6:00 am	88 @ 9:10 am	88 @ 12:53 pm	97 @ 5:53 pm	100 @ 7:24 pm	178 @ 10:50 pm
5/23/2016	4.6 @ 7:30 am	100 @ 8:20 am	100 @ 8:20 am	104 @ 4:36 pm	100 @ 7:10 pm	110 @ 10:00 pm	100 @ 10:00 pm
5/24/2016	4.6 @ 7:30 am	83 @ 7:30 am	79 @ 9:20 am	104 @ 4:36 pm	100 @ 7:10 pm	110 @ 10:00 pm	100 @ 10:00 pm
5/25/2016	4.6 @ 7:30 am	48 @ 7:14 am	44 @ 8:27 am	155 @ 1:14 pm	175 @ 6:57 pm	144 @ 8:11 pm	204 @ 11:23 pm
5/26/2016	4.6 @ 7:30 am	85 @ 8:30 am	104 @ 8:21 am	104 @ 8:21 am	104 @ 8:21 am	104 @ 8:21 am	208 @ 11:23 pm
5/27/2016	4.6 @ 7:30 am	57 @ 8:12 am	100 @ 12:30 pm	100 @ 12:30 pm	79 @ 6:20 pm	102 @ 11:20 pm	102 @ 11:20 pm
5/28/2016	4.6 @ 7:30 am	53 @ 8:00 am	100 @ 12:30 pm	100 @ 12:30 pm	27 @ 6:10 pm	102 @ 11:20 pm	102 @ 11:20 pm
5/29/2016	4.6 @ 7:30 am	53 @ 8:00 am	100 @ 12:30 pm	100 @ 12:30 pm	234 @ 6:50 pm	102 @ 11:20 pm	102 @ 11:20 pm
5/30/2016	4.6 @ 7:30 am	102 @ 10:51 am	111 @ 11:53 am	132 @ 6:12 pm	132 @ 6:12 pm	102 @ 11:20 pm	102 @ 11:20 pm
5/31/2016	4.6 @ 7:30 am	115 @ 8:40 am	115 @ 8:40 am	115 @ 8:40 am	115 @ 8:40 am	346 @ 7:59 pm	140 @ 10:48 pm
6/1/2016	4.6 @ 7:30 am	134 @ 7:00 am	134 @ 7:00 am	134 @ 7:00 am	134 @ 7:00 am	230 @ 7:00 pm	141 @ 11:14 pm
6/2/2016	4.6 @ 7:30 am	140 @ 8:40 am	140 @ 8:40 am	177 @ 1:42 pm	83 @ 4:30 pm	114 @ 10:00 pm	114 @ 10:00 pm
6/3/2016	4.6 @ 7:30 am	127 @ 8:40 am	127 @ 8:40 am	127 @ 8:40 am	70 @ 5:50 pm	70 @ 5:50 pm	101 @ 10:00 pm

Prepare for your next visit with your healthcare professional

- Click **Visiting your doctor?**
The **Be prepared for your visit** screen will appear.
- Click **Download reports**
To include reports, be sure to click the box next to **Include reports for your doctor** before you click **Download reports**.

If you'd like to change the list of reports you want to print for your visit, you can do so by clicking on **Settings**, then **Preferred Reports to Print** and then selecting the reports you need.

Any questions? We're here to help



Visit OneTouch.com



Call us at 1 866 693-0599

**Find the necessary system requirements at
[OneTouchReveal.com/support/en_US/](https://www.onetouchreveal.com/support/en_US/)**

